

# Matrix Retainers

High-Quality Finish in Stainless Steel

**Ergonomic handling**

is achieved through a ball bearing mechanism of the screw, made of PEEK.

**Rounded edges**

minimize the risk of injury.


**The matrix retainer**

securely holds the ends of the matrix bands.




## Nyström

REF 1




## Nystöm-Form

REF 20: 1 Pair  
REF 20XI: I = left  
REF 20XII: II = right




## Ivory-Form

REF 12 with shields



## Ivory-Form

with protection shield  
REF 8X5: Width 5  
REF 8X6: Width 6  
REF 8X7: Width 7



## Tofflemire-Form Universal

REF 21 with removable sliding carriage  
for easy cleaning  
REF 25 with permanent fixed  
sliding carriage

