

Matrixes

Stainless Steel

Micro thin:
0,03 / 0,04 / 0,05 mm

Outstanding tensile strength

HAHNENKRATT matrixes and matrix bands are produced in various shapes. The stainless steel offers outstanding tensile strength.

Practical dispenser

for matrix bands on a roll. The top opening of the dispenser allows for easy sliding out of the matrix band, and if necessary, the band can also be easily pushed back into the dispenser.



Nyström 0,04 mm





Drawings on a scale of 1:1

Matrixes

Ivory 0,04 mm



Tofflemire 0,04 mm



Matrixes



| Matrixes 0,04 mm | REF | | | | | | |
|------------------|--------|--------|--------|--------|--------|--------|------|
| box of 12 | 40X1 | 40X2 | 40X3 | 40X4 | 40X5 | 40X6 | 40X7 |
| | 40X11 | 40X12 | 40X13 | 40X14 | 40X15 | 40X16 | |
| | 40X101 | 40X102 | 40X103 | 40X113 | 40X150 | 40X151 | |
| box of 100 | 41X1 | 41X2 | 41X3 | 41X4 | 41X5 | 41X6 | 41X7 |
| | 41X11 | 41X12 | 41X13 | 41X14 | 41X15 | 41X16 | |
| | 41X101 | 41X102 | 41X103 | 41X113 | 41X150 | 41X151 | |

| Matrix bands in cardboard dispenser | | | | | | | | | | |
|-------------------------------------|--------|--------|-------|-----------------------------|--------|--------|-------|-----------------------------|--------|--------|
| roll of 1m thickness in mm | | | width | roll of 3 m thickness in mm | | | width | roll of 10m thickness in mm | | |
| 0,03 | 0,04 | 0,05 | | 0,03 | 0,04 | 0,05 | | 0,03 | 0,04 | 0,05 |
| 10 134 | 10 144 | 10 154 | 4 mm | 10 334 | 10 344 | 10 354 | 4 mm | 10 034 | 10 044 | 10 054 |
| 10 135 | 10 145 | 10 155 | 5 mm | 10 335 | 10 345 | 10 355 | 5 mm | 10 035 | 10 045 | 10 055 |
| 10 136 | 10 146 | 10 156 | 6 mm | 10 336 | 10 346 | 10 356 | 6 mm | 10 036 | 10 046 | 10 056 |
| 10 137 | 10 147 | 10 157 | 7 mm | 10 337 | 10 347 | 10 357 | 7 mm | 10 037 | 10 047 | 10 057 |

