

# Photography Mirrors

Rhodium Front Surface

ULTRA Front Surface

The precisely sharp reflection is high-contrast down to the finest detail. Combined with brilliant color accuracy and pleasantly bright reflection, you get intraoral images of great clarity.



The capture of the mirror image reveals new perspectives.



**Rhodium FS**  
 Canon 7D and Metz flash:  
 Aperture 25 – Shutter speed  
 125 – ISO 100



**ULTRA FS**  
 The unique construction of the reflective layers – consisting of 22 layer packages – provides an exceptionally bright mirror image through optimized light refraction.



# Photography Mirrors



REF 6660 KO / empty



REF 6530 X ... / 6730 X ...



REF 6540 X ... / 6740 X ...



REF 6550 X ... / 6750 X ...

Form 1, 2, 3, 4 without handle:  
REF 6500 X ... / Rhodium mirrors  
REF 6700 X ... / ULTRA mirrors

Photography Mirrors / each		double-sided without handle	handle bent			handle V-shaped		handle straight	
Form	mm	REF FS Rhodium	REF FS Rhodium	REF ULTRA FS	REF FS Rhodium	REF ULTRA FS	REF FS Rhodium	REF ULTRA FS	
1		58x76	6300 X 1	6530 X 1	6730 X 1			6550 X 1	6750 X 1
		66x95	6300 X 5	6530 X 5	6730 X 5			6550 X 5	6750 X 5
2		74x106	6300 X 2	6530 X 2	6730 X 2			6550 X 2	6750 X 2
3		42x128	6300 X 3	6530 X 3	6730 X 3	6540 X 3	6740 X 3		
4		36x120		6530 X 3KFO	6730 X 3KFO	6540 X 3KFO	6740 X 3KFO		
		55/40x157	6300 X 4						
		53/31x142	6300 X 6						
		86/77x133	6300 X 7						
		71/62x133	6300 X 8						
		47x153	6300 X 9						
Suitcase for Photography Mirrors / empty			REF						
For 3 Photography Mirrors with handle, empty Outside: 285x240x100 mm			6660 KO						

