

# ERGOform

Stainless Steel, satin finished

PEEK

## Ergonomic design

These handles are lightweight and balanced in the hand due to their slim end. The larger diameter and grip recesses allow for a comfortable grip and relaxed working.

## Pleasant haptic

The micro-finely textured, satin-matt surface and gently rounded edges provide a very pleasant feel.

## Satin-finished surfaces


These prevent unwanted light reflections and give the handle a refined finish.

## The ERGOform handle, made of satin-finished stainless steel







is hollow and feels well-balanced and lightweight in the hand.

## The ERGOform handles made of PEEK

PEEK is dimensionally and color-stable up to 210°C and can be sterilized both in an autoclave and with hot air.





| Handles for Mouth Mirrors, Length 135 mm  |             |      | SS-thread<br>REF | CS-thread <sup>*</sup><br>REF |
|---|-------------|------|------------------|-------------------------------|
| ERGOform, PEEK  |             |      | REF              |                               |
|    | ocean blue  | each | 580              | 590                           |
|    | creme white | each | 581              | 591                           |
|    | turquoise   | each | 582              | 592                           |
|    | violet      | each | 584              | 594                           |
|  | anthracite  | each | 585              | 595                           |
| ERGOform, Stainless Steel, hollow   |             |      |                  |                               |
|  |             | each | 492              |                               |

\*Due to a larger thread-diameter, the CS-thread fits only on mouth mirrors with CS-thread.  
CS is mainly used in the USA, Canada, South-Korea and Spain.