

Handles

Stainless Steel polished

Stainless Steel satin finished

The perfect handle for every dentist

Pleasant haptics

The high quality is characterized by the careful finish. The surfaces are finely polished and all edges are smoothly rounded, ensuring the instrument feels comfortable in the hand.

The hollow Handles are well-balanced and lightweight. A wider diameter of some handles and especially the grip recesses allow for a relaxed grip and easy handling.

Satin-finished surface














eliminates disturbing glares and gives the instruments an exclusive finish.

Outstanding stability against corrosion

According to our experience of over 60 years, the stainless steel is stable against corrosion. It complies with EN ISO 7153-1 (Surgical Instruments – Metallic materials).



Handles Stainless Steel

Handles for Mouth Mirrors, Length 130 mm	Features	Content	SS-thread REF	CS-thread [*] REF
	polished, octagonal	each	413	415
	polished, octagonal, serrated	each	420	421
	polished, round, serrated	each	426	427
	polished, round, serrated	each	5012	5013
	polished, with grip recesses	each	430	431
	polished, hexa- gonal, hollow	each	470	471
	satin-finished, hexagonal, hollow	each	476	477
	polished, hexag., hollow 135 mm	each	480	481
	satin-finished, hexagonal, with scale, (not calibrated)	each	474	
	polished, hollow	each	460	461
	satin-finished surface, hollow	each	466	467
	ERGOtouch	each	442	443
	ERGOform, 135 mm	each	492	

*Due to a larger thread-diameter, the CS-thread fits only on mouth mirrors with CS-thread.
CS is mainly used in the USA, Canada, South-Korea and Spain.

