

Economy

The economically priced

High image quality


The Economy mouth mirrors have a pure silver coating. The reflection is bright and offers high image quality.

High resistance against disinfection and sterilisation

The silver coating is protected with a copper layer and it's additionally sealed with a special protective film.



Economy



Economy / Box of 12		Simple Stem (SS)		Cone Socket (CS) *	
Size	Ø in mm	plane	magnifying	plane	magnifying
		REF			
4	22	770X4	771X4	775X4	776X4
5	24	770X5	771X5	775X5	776X5

*Due to a larger thread-diameter, the CS-thread fits only on handles with CS-thread.
CS is mainly used in the USA, Canada, South-Korea and Spain.

