



SEplus

Since decades dentists all over the world rely on the superior quality of our SEplus. Mouth mirrors are daily sterilised and disinfected several times. Considering these extreme conditions, the superior quality is decisive.

Best reflection and image quality


The SEplus coating with a pure silver layer guarantees exceptionally high reflection. The mirror image impresses with distortion-free, color-natural, and bright image quality, enabling fatigue-free viewing in this combination.

Best resistance against disinfection and sterilisation

The silver coating is protected by a combination of layers, made out of copper and titanium, developed by our company. Additionally the mirror is sealed with a special protective film. This unique combination ensures the extreme long life of SE plus mouth mirrors.



SEplus



SEplus / Box of 12		plane	magnifying	plane	magnifying
Size	Ø in mm	Simple Stem (SS)		Cone Socket (CS) *	
		REF			
3	20	702X3	706X3	711X3	716X3
4	22	702X4	706X4	711X4	716X4
5	24	702X5	706X5	711X5	716X5

*Due to a larger thread-diameter, the CS-thread fits only on handles with CS-thread.
CS is mainly used in the USA, Canada, South-Korea and Spain.

