

MEGA FS Rhodium

Form follows function –
it's ergonomically designed

More comfortable cheek retraction
due to the flat, broad-shaped stem.


No disturbing reflection
by satin-finished stainless steel surface.

Exact sharpness
by Front Surface mirror coating

Easy cleaning
due to the anti-stick effect of the special
HAHNENKRATT front surface coating.



MEGA FS Rhodium



MEGA FS Rhodium				SS-thread REF	CS-thread REF	
1	MEGA	4 plane	Ø 22 mm	box of 12	6 060	6 070
	MEGA	5 plane	Ø 24 mm	box of 12	6 061	6 071
2	MEGAduo, double sided	4 plane	Ø 22 mm	box of 6	6 100	6 170
	MEGAduo, double sided	5 plane	Ø 24 mm	box of 6	6 101	6 171
3	MEGAretract	4 plane	Ø 22 mm	box of 6	6 660	6 665
	MEGAretract	5 plane	Ø 24 mm	box of 6	6 661	6 666
4	MEGAretract, bifurcated stem	4 plane	Ø 22 mm	box of 6	6 680	6 685
	MEGAretract, bifurcated stem	5 plane	Ø 24 mm	box of 6	6 681	6 686
5	MEGAstretch		18x34 mm	each	6 860	
	MEGAstretch		22x38 mm	each	6 861	
6	MEGA Parallelometer		22x37 mm	each	6 900	
	MEGA Wangenabhalter	without mirror	Ø 22 mm	each	6 070 W	
	MEGAabrasiv	Saphire-hard glass surface suitable for use with aluminium oxide airbrasive	4 plane	Ø 21 mm	each	6 080

*Due to a larger thread-diameter, the CS-thread fits only on handles with CS-thread.
CS is mainly used in the USA, Canada, South-Korea and Spain.

