


ULTRAvision FS

Front Surface


The unique structure of the reflective layers – consisting of 22 layer packages – offers an unmatched bright reflection with crystal-clear purity through the optimized refraction of light, which is not found in other mirrors.

Highest reflection
with crystal clear brightness




Superior colour accuracy,

Exact sharpness
by Front Surface mirror coating



ULTRAvision FS



ULTRAvision FS			Simple Stem(SS) plane	Cone Socket (CS) plane
Size	Ø in mm		REF	REF
0	14	box of 3	532X0	
3	20	box of 12	533X3	536X3
4	22	box of 12	533X4	536X4
5	24	box of 12	533X5	536X5

*Due to a larger thread-diameter, the CS-thread fits only on handles with CS-thread.
CS is mainly used in the USA, Canada, South-Korea and Spain.