

ULTRA FS

ultra bright

ULTRA FS – the brightest and sharpest mirror worldwide.

Highest reflection
with crystal clear brightness

Superior colour accuracy

Images of exact sharpness and contrast
up to the smallest detail, with out any distortion.
The combination of bright and brilliant natural colours guarantees high sight comfort, gentle to the eye and good for fatigue-proof working.

Optimized comfortable cheek retraction
through wide, contoured stem bases




ULTRA FS

ultra bright

The unique structure of the reflective layers – consisting of 22 layer packages – offers an unmatched bright reflection with crystal-clear purity through the optimized refraction of light, which is not found in other mirrors.


only one reflection
▶ exact sharp image



Front Surface Mirror



The splitted mirror shows the crystal clear, ultra-bright reflection of the ULTRA FS compared to Rhodium.



ULTRA FS



				SS-thread REF	CS-thread [*] REF
1	ULTRA FS	4 plane	box of 6	6 160	6 165
	ULTRA FS	5 plane	box of 6	6 161	6 166
2	ULTRAduo FS, double sided	4 plane	box of 6	6 640	6 645
	ULTRAduo FS, double sided	5 plane	box of 6	6 641	6 646
3	ULTRAretract FS	4 plane	box of 6	6 670	6 675
	ULTRAretract FS	5 plane	box of 6	6 671	6 676
4	ULTRAretract FS, bifurcated stem	4 plane	box of 6	6 690	6 695
	ULTRAretract FS, bifurcated stem	5 plane	box of 6	6 691	6 696
5	ULTRAretract FS, bifurcated stem stem double 45°	4 plane	each	6 650	

*Due to a larger thread-diameter, the CS-thread fits only on handles with CS-thread.
CS is mainly used in the USA, Canada, South-Korea and Spain.