

# BLACKline

The BLACKline instruments receive with a DLC (Diamondlike carbon) coating a high-end finish.

The matt finish eliminates disturbing reflections. This is an important aspect for working under the microscope, which is becoming increasingly important in endodontics.

The BLACKline increases working comfort while adding ergonomic elements and offers superior advantages for every purpose.

# BLACKline

## Handles + Explorers


The handles with DLC coating are hollow, making them perfectly balanced and lightweight in the hand. The textured design and recess provide a secure grip and pleasant tactile experience. The deep black, matt surface eliminates distracting light reflections.


The explorer parts are made out of surgical steel according to EN ISO 7492 and thus ensure a high bending strength and flexibility.







The handles are also available for mouth mirrors.

**The Explorers are permanently fixed.**


**Autoclavable up to 134°C**





BLACKline ERGOtouch Explorers / each	REF				
	 <b>Form A</b> 010 Fig. 23	 <b>Form C</b> 030 Fig. 8	 <b>Form D</b> 040 Fig. 9	 <b>Form F</b> 060 Fig. 6	 <b>Form G</b> 070 Fig. 17
	B600A	B600C	B600D	B600F	B600G

The Explorers are permanently fixed.

BLACKline ERGOtouch Handles / each	SS-thread	CS-thread*
	B442	B443

\*Due to a larger thread-diameter, the CS-thread fits only on handles with CS-thread. CS is mainly used in the USA, Canada, South-Korea and Spain.