



Years of  
Passion for Dental  
Excellence



Dental Excellence  
made by

**HAHNENKRATT**

Made in Königsbach-Stein



1932

Emil  
Hahnenkratt

## It began more than 90 years ago.

The company E. HAHNENKRATT GmbH was founded in 1932 by master mechanic Emil Hahnenkratt in Pforzheim.

The passion for technology has been passed down through all generations, even to the great-grandchild, and is now reflected in the consistent use of high-tech in all areas of the company.



Ann-Katrin and Hendrik Schepperheyne

2025

## WE LOVE TO WELCOME YOU.

Say HELLO to our IDS NEWcomers 2025:

**CONTECpro**

**CYTECpro**

**LightTransmitting**

It would be great if we could excite you for our new post systems and we are looking forward to hear your feedback.

Your **HAHNENKRATT Team**  
with the Schepperheyne family

Ann-Katrin Schepperheyne

– CSO –

[a.schepperheyne@hahnenkratt.com](mailto:a.schepperheyne@hahnenkratt.com)



For information and inspiration, please visit  
[www.hahnenkratt.com](http://www.hahnenkratt.com)





## EXATEC CONTEC CYTEC

Scientifically tested  
and proven in practice  
since 2001.

### High bending strength

Achieved through the use of high tenacity HT Glassfiber, which bonds perfectly with the matrix:

HT Glassfiber: 1.678 MPa  
Determined accordingly to EN ISO 178

### High fatigue and fracture resistance

Homogeneity due to dentin-like elastic modulus

Secure adhesion due to micro-retentive surface

High radiopacity

## CONTEC<sup>pro</sup> CYTEC<sup>pro</sup>

### LightTransmitting

The new **pro** root posts are a good addition to the EXATEC, CYTEC, and CONTEC root posts made of white opaque HT-Glassfiber, which have been proven in practice worldwide for over 20 years.

The new **pro line** features the same excellent properties, particularly in terms of stability and secure bonding with composite materials.

**pro** LightTransmitting throughout the entire apical post area for – optimized and shortened – primary curing of the composite, up to the tip of the root post. In a tightly timed schedule, the dentist gains valuable seconds to be able to proceed with the next step without losing time

Get **YOUR**  
free IDS promo set

[aftersales@hahnenkratt.com](mailto:aftersales@hahnenkratt.com)

IDS 2025  
NEWS

## BLACKline

ULTRA FS  
MEGA FS Rhodium

»Diamond Like Carbon« DLC coating  
for High-End Finish.

**Eye-friendly work**  
under bright light and magnification due to the  
deep black, matt surface.

**Precisely sharp reflection**  
thanks to front surface coating.

**Highest reflection**  
for brilliant brightness, especially with ULTRA FS.

**Optimized cheek retraction**  
through wide, contoured stem bases.



## RELAX

created with a perfect blend  
of design and ergonomics.

**Precisely sharp reflection**  
thanks to front surface coating.

**Highest reflection**  
for brilliant brightness, especially with ULTRA FS.

**Comfortably light weight**  
for relaxed dental operations

**Comfortable to hold**  
due to the ergonomic design with grip recess

**Keeping quality**  
colour stable, resistant against acids as well  
as against plaque-indicators

**Hygienic**  
due to the plane transition from mirror surface  
to mirror casing without any gap

**No distracting reflections**  
due to the satin-finished surfaces

**The instrument is autoclavable, acid-resistant  
and suitable for thermal disinfectors.**





**E. HAHNENKRATT GMBH**  
**Dentale Medizintechnik**

---

Benzstraße 19  
DE-75203 Königsbach-Stein

Fon +49 7232 3029-0

info@hahnenkratt.com  
sales@hahnenkratt.com  
marketing@hahnenkratt.com

www.hahnenkratt.com

